

idi**Aspect High-Tilt™ Mobile Imaging Table**

Product Data Sheet

*Patent pending***Ideal C-arm imaging solution for:**

- ***Swallowing studies for post-surgical gastric banding adjustment***
- ***Other C-arm procedures requiring upright patient imaging or loading***

- Motorized tilt from horizontal to 75° upright with footswitch operator control
- “Step-on” patient loading with table in upright position
- 525 lb. patient weight capacity
 - Ultra-stable during tilt motion and in any table position
- Easy C-arm imaging access in any tabletop position
 - Compact table pedestal = very large imaging area
- Carbon fiber tabletop: artifact-free, low x-ray attenuation
- Standard 110 volt operation with battery back-up
- 5”, locking casters and lightweight table for easy transport
- Table includes:
 - built-in footboard
 - tabletop pad
 - patient handgrips
 - restraint straps
 - footswitch control



Aspect High-Tilt™ Product Specifications

- Table tilt range: 0° (horizontal) to 75° upright
- Tilt speed: $\approx 4^\circ$ per second, ≈ 20 seconds to full stroke
- Tabletop height in horizontal position: 32" (81cm)
- Footboard position with table upright: 5" (12.7cm) above floor
- Tabletop size: 24"W x 80"L (61cm x 203cm)
- Total imaging area: ≈ 55 "L (140 cm) x 24" W (61cm)
- Standard US siderails ($\frac{3}{8}$ " x $1\frac{1}{8}$ "") for handgrips attachment
- Patient weight capacity: 525 lb. (238 kg)
- Table weight: ≈ 425 lbs. (192.8kg)
- Electrical power: 110V / ≤ 10 amp

