



## Aspect 100US Gyn/Fertility Clinic Mobile Imaging Table

*Ideal C-arm Imaging Solution for Hysterosalpingograms*



- Carbon fiber tabletop with low x-ray attenuation for superior fluoroscopic imaging
- Unobstructed fluoroscopy imaging area: 22"W x 32"L (56 x 81 cm)
- Motorized tabletop elevation with footswitch control: 26" – 38" (66 – 96.5cm)
- 24"W x 46"L tabletop (61 x 117 cm)
- Patient weight capacity: 450 lbs. (204 kg)
- Standard US table side-rail: compatible with many lithotomy legholders and clamps.
- 120 VAC operation with battery back-up operation (230VAC optional)
- Easy table transport with 5" locking casters. Table weight = 375 lbs. (170 kg)
- 3-year parts warranty
- Accessories included:
  - tabletop pad, patient restraint straps, footswitch control
- Optional accessories available:
  - Lithotomy legholders and adjustable clamps (pictured above)
  - Table hand control