




Tabletops for Aspect Mobile Imaging Tables


All Aspect Series Tabletops feature very low attenuation, carbon fiber construction.
Optional tabletop styles below must be selected at time of table order.

Aspect 100, 100E, 100T, 100R, 100RT and 100RTL Tables



Standard	Optional
X100-995	X100-733
24" x 80" 61 x 203cm	22" x 80" 56 x 203cm with facial cutout

Aspect 100-4 Table



Standard	Optional	Optional
X100-517	X100-982	X100-1597
18" - 30" W x 80" L 45.7 - 76.2 x 203cm	23.5" - 30" W x 80" L 59.4 - 76.2 x 203cm	18" - 27" W x 80" L 45.7 - 68.6 x 203cm