



**Aspect Series Tables Specifications & Accessory Matrix:  
Pain Management, Vascular and General Imaging Tables**

	100	100E	100R	100T	100RT	100RTL	100-4	100-4T	100EP
Elevate (tabletop to floor)	Fixed	26 - 37.8" 66 - 96 cm	28.5 - 40.3" 72 - 102 cm	28 - 39.8" 71 - 101 cm	30.4 - 42.2" 77 - 107 cm	31.2 - 41.2" 79 - 105 cm	27.9 - 37.7" 71 - 96 cm	29 - 38.9" 73.7 - 99 cm	29 - 38.9" 74 - 99 cm
Trendelenburg Tilt	NA	NA	NA	±15°	±15°	±15°	NA	±10°	75 / 15°
Lateral Tilt	NA	NA	20°	NA	20°	20°	NA	NA	NA
Longitudinal Tabletop Motion	NA	NA	NA	NA	NA	10" 25 cm	40" float 102 cm	32" float 81 cm	NA
Transverse Tabletop Motion	NA	NA	NA	NA	NA	NA	10" float 25 cm	8" float 20 cm	NA
Standard Tabletop Size	24" x 80" 61 x 203 cm	24" x 80" 61 x 203 cm	24" x 80" 61 x 203 cm	24" x 80" 61 x 203 cm	24" x 80" 61 x 203 cm	24" x 80" 61 x 203 cm	18 -30" x 80" 46-76 x 203 cm	24" x 84" 61 x 213 cm	22" x 84" 56 x 213 cm
Optional Tabletop Size	22" x 80" 56 x 203 cm	22" x 80" 56 x 203 cm	22" x 80" 56 x 203 cm	22" x 80" 56 x 203 cm	22" x 80" 56 x 203 cm	22" x 80" 56 x 203 cm	25-30" x 80" 64-76 x 203 cm	NA	NA
<i>Note: Optional 22" x 80" tabletops have bottleneck shape (narrow section at head-end)</i>									
Imaging Area	66" 167.6 cm	66" 167.6 cm	62" 157.5 cm	62" 157.5 cm	62" 157.5 cm	59" 149.9 cm	72" 182.9 cm	62" 157.5 cm	61" 154.9 cm
Patient Capacity	450 lb. 204 kg	450 lb. 204 kg	450 lb. 204 kg	450 lb. 204 kg	450 lb. 204 kg	450 lb. 204 kg	450 lb. 204 kg	450 lb. 204 kg	350 lb. 159 kg
Table Weight	250 lb. 113 kg	290 lb. 132 kg	325 lb. 147 kg	325 lb. 147 kg	325 lb. 147 kg	355 lb. 161 kg	385 lb. 175 kg	500 lb. 227 kg	425 lb. 193 kg
<b>Standard (●) and Optional (O) Accessories</b>									
Tabletop Pad	●	●	●	●	●	●	●	●	●
Hand Control	NA	●	●	●	●	●	●	●	O
Panning Hand Control	NA	NA	NA	NA	NA	NA	●	●	NA
Foot Control	NA	O	O	O	O	O	O	O	●
ID410 Width Extender	Optional. For use with 22"wide tabletops only.						●	NA	NA
Facial Support Pad	O	●	●	●	●	●	NA	NA	NA
Armboard, rail mounted	O	O	O	O	O	O	NA	NA	NA
Armboard, carbon fiber	O	O	O	O	O	O	O	O	O
Restraint Straps	●	●	●	●	●	●	●	●	●
ID400 Dual Armboard	Optional. Compatible with bottleneck shape tabletops only.							O	O
Clamp-on Accessory Rail	O	O	O	O	O	O	NA	NA	NA
Catheter Tray extension	NA	NA	NA	NA	NA	NA	O	NA	NA
Footboard	NA	NA	NA	NA	NA	NA	NA	NA	●
Patient Handgrips	NA	NA	NA	NA	NA	NA	NA	NA	●
Table-side Radiation Shield	O	O	O	O	O	O	O	O	NA



**Aspect Urology Table Specifications & Accessory Matrix**

	<b>100US</b>	<b>100U</b>	<b>100UC</b>
Elevate (tabletop to floor)	26 - 38" ( 66 – 96.5 cm )	31 - 40.9" ( 78.8 – 103.8 cm )	31.27 - 41.1" ( 79.4 – 104.4 cm )
Trendelenburg Tilt	NA	±15°	±15°
Longitudinal Tabletop Motion	NA	10" ( 25.4 cm )	10" ( 25.4cm )
Transverse Tabletop Motion	NA	NA	± 3.5" ( ±8.9cm )
Tabletop Size	24"W x 46"L ( 61 x 117 cm )	28"W x 46"L ( 71 x 117 cm )	28"W x 46"L ( 71 x 117 cm )
Imaging Area	22"W x 32"L ( 56 x 81 cm )	25"W x 19"L (63.5 x 48.3 cm)	25"W x 19"L (63.5 x 48.3cm )
Patient Capacity	450 lb. ( 204 kg )	525 lb. ( 238 kg )	525 lb. ( 238 kg )
Table Weight	≈ 375 lb. ( ≈170 kg )	≈ 575 lb. ( ≈261 kg )	≈ 575 lb. ( ≈261 kg )
<b>Standard (●) and Optional (O) Accessories:</b>			
Tabletop Pad	●	●	●
Hand Control	O	●	●
Foot Control	●	●	●
Armboard, rail mounted	O	●	●
Lithotomy Legholders (pair)	O	O	O
X-ray Cassette Tray	NA	●	●
Uro Drain Bag Hoop	NA	●	●
Urologist Elbow Supports (pair)	NA	●	●
Tabletop Foot Extension	NA	●	●
X-ray shield drape for perineal end of table	NA	●	●

Note: Lithotomy Legholders are optional on all uro table models. See Product Catalog for Options.

**Aspect Series Tables Electrical Specifications & Standards Compliance**

All Aspect Series tables operate on 110 – 120 VAC @ less than 10 A. 230VAC operation is optional and must be specified with order.

All IDI Aspect Series Tables are ETL Listed (IEC, UL & CSA Standards) and CE marked.



Technical Specifications subject to change without notice.

**Image Diagnostics, Inc., 98 Pratt's Junction Rd.. P.O. Box 576, Sterling, MA 01564-0576  
Tel : 978.422.8601**

D000-499G  
4.1.08