



Tables Specifications & Accessory Matrix

Pain Management and General Imaging Tables

	100E	EZ Access	100R	100T	100RT	100RTL	High-Tilt
Elevate*	26 - 37.8" 66 - 96 cm	21 - 40" 53 - 102 cm	28.5 - 40.3" 72 - 102 cm	28 - 39.8" 71 - 101 cm	30.4 - 42.2" 77 - 107 cm	31.2 - 41.2" 79 - 105 cm	Fixed: 32" 81 cm
Trendelenburg Tilt	NA	NA	NA	±15°	±15°	±15°	75° / 0°
Lateral Tilt	NA	NA	20°	NA	20°	20°	NA
Longitudinal Tabletop Motion	NA	NA	NA	NA	NA	10" 25 cm	NA
Transverse Tabletop Motion	NA	NA	NA	NA	NA	NA	NA
Standard Tabletop Size	24" x 80" 61 x 203 cm	24" x 80" 61 x 203 cm	24" x 80" 61 x 203 cm	24" x 80" 61 x 203 cm	24" x 80" 61 x 203 cm	24" x 80" 61 x 203 cm	24" x 80" 61 x 203 cm
Optional Tabletop Size	22" x 80" 56 x 203 cm	22" x 80" 56 x 203 cm	22" x 80" 56 x 203 cm	22" x 80" 56 x 203 cm	22" x 80" 56 x 203 cm	22" x 80" 56 x 203 cm	NA
Imaging Area	66" 167.6 cm	68" 172.7 cm	62" 157.5 cm	62" 157.5 cm	62" 157.5 cm	59" 149.9 cm	61" 154.9 cm
Patient Capacity	450 lb. 204 kg	450 lb. 204 kg	450 lb. 204 kg	450 lb. 204 kg	450 lb. 204 kg	450 lb. 204 kg	525 lb. 238 kg
Table Weight	350 lb. 159 kg	450 lbs. 208 kg	395 lbs. 179 kg	395 lbs. 179 kg	425 lbs. 193 kg	460 lb. 209 kg	605 lb. 274 kg
Tabletop Pad	●	●	●	●	●	●	●
Restraint Straps	●	●	●	●	●	●	●
Facial Support Cushion	●	●	●	●	●	●	NA
Hand Control	●	●	●	●	●	●	NA
Foot Control	○	NA	○	○	○	○	●
ID410 Width Extender	○ For 22" top only	○ For 22" top only	○ For 22" top only	○ For 22" top only	○ For 22" top only	○ For 22" top only	NA
Armboard, rail mounted	○	○	○	○	○	○	NA
Armboard, carbon fiber	○	○	○	○	○	○	○
Clamp-on Accessory Rail	○	○	○	○	○	○	NA
Patient Handgrips	NA	NA	NA	NA	NA	NA	●
Table Mounted Radiation Shield	○	○	○	○	○	○	NA

* Table elevation is distance from tabletop to floor *not* including 2" (2.5cm) tabletop pad



Tables Specifications & Accessory Matrix

Tables with Manual 4-way Floating Tabletops

	100-4	100-4T	ISR
Elevate*	27.9 - 37.7" 71 - 96 cm	29 - 38.9" 73.7 - 99 cm	30 - 38" 76 - 96.5 cm
Trendelenburg Tilt	NA	±10°	±12°
Lateral Tilt	NA	NA	±12°
Longitudinal Tabletop Motion	40" float 102 cm	32" float 81 cm	32" float 81 cm
Transverse Tabletop Motion	10" float 25 cm	8" float 20 cm	8" float 20 cm
Standard Tabletop Size	18 -30" x 80" 46-76 x 203 cm	24" x 84" 61 x 213 cm	22" - 24" x 84" 56-61 x 213 cm
Optional Tabletop Size	23.5-30" x 80" 59-76 x 203cm	NA	24" x 84 61 x 213 cm
Imaging Area	72" 182.9 cm	62" 157.5 cm	61", 73" with optional head extension 155 cm, 185 cm with optional head extension
Patient Capacity	450 lb. 204 kg	450 lb. 204 kg	450 lb. 204 kg
Table Weight	490 lb. 222 kg	605 lb. 274 kg	605 lb. 274 kg
Tabletop Pad	•	•	•
Restraint Straps	•	•	•
Panning Hand Control, rail mounted	•	•	•
Pendant Hand Control	NA	NA	○
Foot Control	○	○	○
Armboard, rail mounted	○	○	○
Armboard, carbon fiber	○	○	○
VAB armboard	○	○	○
Clamp-on Accessory Rail	○	○	○
Catheter Tray extension	○	○	○
Head Extension, Articulating	NA	NA	○
Anesthesia Screen Holder	○	○	○
Table Mounted Radiation Shield	○	○	○

* Table elevation is distance from tabletop to floor *not* including 2" (2.5cm) tabletop pad



**Tables Specifications & Accessory Matrix
Urology and GYN Tables**

	100US	100UC Plus
Elevate (tabletop height, excluding tabletop pad thickness)	26 - 38" (66 – 96.5 cm)	31.27 - 41.1" (79.4 – 104.4 cm)
Trendelenburg Tilt	NA	±15°
Longitudinal Tabletop Motion	NA	10" (25.4cm)
Transverse Tabletop Motion	NA	± 3.5" (±8.9cm)
Tabletop Size	24"W x 46"L (61 x 117 cm)	28"W x 46"L (71 x 117 cm)
Uro/Gyn Imaging Area	22"W x 32"L (56 x 81 cm)	25"W x 19"L (63.5 x 48.3cm)
Tabletop Leg Extension	NA	24"W x 34"L (71 x 86.4 cm)
Total Imaging Area with Extension	NA	24"W x 53"L (61 x 135 cm)
Patient Capacity	450 lb. (204 kg)	525 lb. (238 kg)
Table Weight	≈ 395 lb. (≈179 kg)	≈ 670 lb. (≈304 kg)
Tabletop Pad	●	●
Hand Control	○	●
Foot Control	●	●
Armboard, rail mounted	○	●
Lithotomy Legholders (pair)	○	○
X-ray Cassette Tray	NA	●
Uro Drain Bag Hoop	NA	●
Urologist Elbow Supports	NA	○
X-ray shield drape for perineal end of table	NA	●

Note: Lithotomy Legholders are optional on all uro table models. See Product Catalog for Options.

Aspect Series Tables Electrical Specifications

All Aspect Series tables operate on 110 – 120 VAC @ less than 10 Amps. 230VAC operation is optional and must be specified with order.

Technical Specifications subject to change.

Image Diagnostics, Inc., Fitchburg, MA, USA

www.imagediagnostics.com